

ESA612

Electrical Safety Analyzer

Technical Data



The ESA612 Electrical Safety Analyzer represents the next generation in testers for biomedical professionals that perform field service on medical equipment throughout their facilities, in clinics, and anywhere onsite service is required. Portable, lightweight, and designed for operation in tight spaces, the ESA612 offers the functionality of a simulator, multimeter and electrical-safety analyzer in a single test tool.

With selection of two test loads, this versatile product can be used worldwide to test to preventative maintenance electrical safety standards of choice: ANSI/AAMI ES1:1993 (NFPA-99), IEC62353 (VDE 751), and AN/NZS 3551.

The versatility of the multifaceted ESA612 is further expanded with optional automation software, which speeds and simplifies testing and provides high-end-analyzer productivity at software-level investment. Ansur-automated ESA612 standardizes test procedures, compares results to standards limits, and generates and stores reports for total digital data management.

Key features

- Portable, ergonomic, lightweight and easy to use
- Large, easy-to-read display with adjustable contrast
- Human-factors-designed user interface
- Tilt stand design for stand-up testing in field environments
- Five applied parts jacks and easy ECG snap connection with optional expander box
- ECG waveform tests and dual-lead measurements combine the functionality of a simulator, multimeter and electrical-safety analyzer in a single test tool
- Replaceable mains fuses keep the device in the field and out of the repair shop
- Internal memory for 100 test records
- 20 A at 120 V current capability
- USB connection for use with Ansur and Data Viewer software (for memory download to PC)
- Two-year extended warranty (no-cost, available after first-year calibration at the Fluke Biomedical Cleveland Service Center)
- Optional Ansur automation software standardizes test procedures, compares results to standards limits, generates/stores reports and provides total digital data management
- Rigorously tested for rugged field applications, with CE and CSA in addition to the Fluke-quality-design stamp of approval

Specifications

| Voltage | | |
|---|--|--|
| Range (mains voltage) | 90 V ac to 132 V ac rms | |
| | 180 V ac to 264 V ac rms | |
| Range (accessible voltage) | 0 V ac to 300 V ac rms | |
| Accuracy | ± (2 % of reading + 0.2 V) | |
| Voltage tests | Mains and point-to-point | |
| Earth resistance | | |
| Mode | Two terminal | |
| Test current | > 200 mA ac | |
| Range | 0 Ω to 2 Ω | |
| Accuracy | ± (2 % of reading + 0.015 Ω) | |
| Resistance tests | Earth resistance and point-to-point | |
| Equipment current | | |
| Mode | AC rms | |
| Range | 0 A to 20 A | |
| Accuracy | ± 5 % of reading + (2 counts or 0.2 A, whichever is greater) | |
| Duty cycle | 15 A to 20 A, 5 min on/5 min off 10 A to 15 A, 7 min on/3 min off 0 A to 10 A continuous | |
| Leakage current | | |
| Modes* | AC + DC (true-rms) | |
| | AC only | |
| | DC only | |
| *Modes are available in all leakage tests with the exception of MAP leakages that are available only in true-rms. | | |
| Patient load selection (input impedance) | AAMI ES1-1993 Fig.1 | |
| | IEC 60601-1: Fig 15 | |
| Crest factor | ≤ 3 | |
| Ranges | 0 μA to 199.9 μA | |
| | 200 μA to 1999 μA | |
| | 2 mA to 10 mA | |
| Frequency response/accuracy | DC to 1 kHz | ± (1 % of reading + (1 μA or 1 LSD, whichever is greater)) |
| | 1 kHz to 100 kHz | ± (2 % of reading + (1 μA or 1 LSD, whichever is greater)) |
| | 1 kHz to 5 kHz (current > 1.6 mA) | ± (4 % of reading + (1 μA or 1 LSD, whichever is greater)) |
| | 100 kHz to 1 MHz | ± (5 % of reading + (1 μA or 1 LSD, whichever is greater)) |
| Note: Accuracy for Isolation, MAP, Direct AP, Alternative AP, and Alternative Equipment leakage tests all ranges are + (2.5 μA or 1 LSD, whichever is greater) | | |
| Leakage tests | Ground wire (earth) | |
| | Chassis (enclosure) | |
| | Lead to ground (patient) | |
| | Lead to lead (patient auxiliary) | |
| | Lead isolation (mains on applied part) | |
| | Direct equipment | |
| | Direct applied part | |
| | Alternative equipment | |
| | Alternative applied part | |
| | Point to point | |

| | | |
|---|---|---------------------------|
| Mains on applied part test voltage | 100 % of mains | |
| Differential leakage | | |
| Ranges | 75 μ A to 199 μ A | |
| | 200 μ A to 2000 μ A | |
| | 2 mA to 20 mA | |
| Accuracy | \pm 10 % of reading + (2 counts or 20 μ A, whichever is greater) | |
| Insulation resistance | | |
| Ranges | 0.5 M Ω to 20 M Ω | |
| | 20 M Ω to 100 M Ω | |
| Accuracy | \pm (2 % of reading + 0.2 M Ω) | |
| | \pm (7.5 % of reading + 0.2 M Ω) | |
| Source test voltage | 500 V dc | |
| | 250 V dc | |
| Insulation resistance tests | Mains-PE, AP-PE, Mains- PE, Mains-NE (non-earthed accessible conductive part) and AP- NE (non-earthed accessible conductive part) | |
| ECG performance waveforms | | |
| Accuracy | \pm 2 % | |
| | \pm 5 % for amplitude of 2 Hz square wave only, fixed at 1 mV Lead II configuration | |
| Waveforms | Rates | |
| | ECG complex (BPM) | 30, 60, 120, 180, and 240 |
| | Ventricular fibrillation | |
| | Square wave (50 % duty cycle) (Hz) | 0.125 and 2 |
| | Sine wave (Hz) | 10, 40, 50, 60, and 100 |
| | Triangle wave (Hz) | 2 |
| Pulse (63 ms pulse width) | 30 BPM and 60 BPM | |
| Power ratings | | |
| Mains voltage outlet | 120 V ac or 230 V ac | |
| Mains voltage inlet power range | 90 to 132 V ac rms | 180 to 264 V ac rms |
| Maximum current | 20 A | 16 A |
| Hz | 50 or 60 | 50 or 60 |
| Physical case | | |
| Dimensions (L x W x H) | 17.63 cm x 8.38 cm x 28.45 cm (6.94 in x 3.30 in x 11.20 in) | |
| Weight | 1.6 kg (3.5 lb) | |
| Environmental specifications | | |
| Operating temperature | 10 $^{\circ}$ C to 40 $^{\circ}$ C (50 $^{\circ}$ F to 104 $^{\circ}$ F) | |
| Storage temperature | -20 $^{\circ}$ C to 60 $^{\circ}$ C (-4 $^{\circ}$ F to 140 $^{\circ}$ F) | |
| Operating humidity | 10 % to 90 % non-condensing | |
| Altitude | 120 V ac mains supply voltage up to 5,000 m | |
| | 230 V ac mains supply voltage up to 2,000 m | |
| General | | |
| Warranty | Two-year extended warranty (no-cost, available after first-year calibration at the Fluke Biomedical Cleveland Service Center, otherwise standard one-year warranty applies) | |

Ordering information

Item numbers/Descriptions

ESA612 Electrical Safety Analyzer

ESA612 ESA612 United States, 115 V, 20 A

ESA612-01 ESA612 France, 230 V

ESA612-02 ESA612 Europe, 230 V

ESA612-03 ESA612 Israel, 230 V

ESA612-05 ESA612 Australia, 230 V

ESA612-06 ESA612 United Kingdom, 230 V

ESA612-07 ESA612 Switzerland, 230 V

ESA612-08 ESA612 Thailand, 230 V

ESA612-09 ESA612 Japan, 100 V

ESA612-10 ESA612 North America, 220 V

TA-ESA612-USA ESA612 United States, 115 V, 20 A

with test automation

TA-ESA612-EUR ESA612 Europe, 230 V with test automation

TA-ESA612-FR ESA612 France, 230 V with test automation

TA-ESA612-ISR ESA612 Israel, 230 V with test automation

TA-ESA612-AUS ESA612 Australia, 230 V

with test automation

TA-ESA612-UK ESA612 United Kingdom, 230 V with test automation

TA-ESA612-SWI ESA612 Switzerland, 230 V with test automation

TA-ESA612-THAI ESA612 Thailand, 230 V with test automation

TA-ESA612-JAPAN ESA612 Japan, 100 V with test automation

TA-ESA612-NA220V ESA612 North America, 220 V with test automation

Standard accessories

CD-ROM Operator's Manual (multilingual CD)

MANUAL Getting-Started Guide (hard copy, multilingual)

CD-ROM Ansur ESA612 Plug-In, CD with demo version

CABLE ASSEMBLY Data Transfer Cable

ESA T/L KIT USA US Accessory Kit (included for US, Australia, Israel, Thailand, and Japan versions only):

– Test Lead Set

– TP1 Test Probe Set

– AC285 Alligator Clip Set

ESA T/L KIT EUR EUR Accessory Kit (included for Europe, France, United Kingdom, and Switzerland versions only)

– Test Lead Set

– TP74 Test Probe Set

– AC285 Alligator Clip Set

2719-0154 15 A to 20 A Adapter (US only)

ESA620-NPA Null Post Adapter

ESA612-2016 5-to-5 Banana Jack to ECG (BJ2ECG) Adapter

9530-0075 Carry Case

Power Cord, one included, country specific by model number

LINE CORD US

75026 Europe

75024 UK

75025 Australia

LINE CORD Israel

75026 France/Belgium

75033 Thailand

75058 Japan

FBC-ESA620-4420 Switzerland

Optional accessories

6358 Retractable Test Leads

9503-0004 Ground Pin Adapter (US receptacle testing ground lug)

1210 ECG 1210 Adapter Box Assembly

ANSUR ESA612 Ansur ESA612 Plug-In License Key

About Fluke Biomedical

Fluke Biomedical is the world's leading manufacturer of quality biomedical test and simulation products. In addition, Fluke Biomedical provides the latest medical imaging and oncology quality-assurance solutions for regulatory compliance. Highly credentialed and equipped with a NVLAP Lab Code 200566-6 accredited laboratory, Fluke Biomedical also offers the best in quality and customer service for all your equipment calibration needs.

Today, biomedical personnel must meet the increasing regulatory pressures, higher quality standards, and rapid technological growth, while performing their work faster and more efficiently than ever. Fluke Biomedical provides a diverse range of software and hardware tools to meet today's challenges.

Fluke Biomedical Regulatory Commitment

As a medical test device manufacturer, we recognize and follow certain quality standards and certifications when developing our products. We are ISO 9001 certified and our products are:

- CE Certified, where required
- NIST Traceable and Calibrated
- UL, CSA, ETL Certified, where required
- NRC Compliant, where required

Fluke Biomedical.

Better products. More choices. One company.

Fluke Biomedical

6045 Cochran Road
Cleveland, OH 44139-3303 U.S.A.

Fluke Biomedical Europe

Science Park Eindhoven 5110
5692EC Son, The Netherlands

For more information, contact us:

In the U.S.A. (800) 850-4608 or
Fax (440) 349-2307

In Europe/M-East/Africa +31 40 267 5435 or
Fax +31 40 267 5436

From other countries +1 (440) 248-9300 or
Fax +1 (440) 349-2307

Email: sales@flukebiomedical.com

Web access: www.flukebiomedical.com

©2009-2010 Fluke Biomedical. Specifications subject to change without notice. Printed in U.S.A.
4/2010 3456744C D-EN-N

Modification of this document is not permitted without written permission from Fluke Corporation.